



EZ Way, Inc.

*"Your Total Patient
Lift Solution"*

Rise & Go™ | Walker with Power Rise

A Cost-effective Walker Enhanced with a patented Power Rise Function and Articulating Legs.

BENEFITS:

- **Cost Efficient** - Rise & Go is a complete solution that saves time, space and money.
- **Safe, Confident** - The Rise & Go provides safe and confident standing and mobility. This creates a better working environment for both patients and health care personnel.
- **Power Rise Assistance** - The unique power rise assistance provides a natural standing action.
- **Washable Harness** - The harness is washable, simple, and flexible to wear even while toileting.
- **Height Adjustment** - Rise & Go is equipped with a smooth electric height adjustment.

FEATURES:

- **Power Rise** for assisting patients while standing.
- **Articulating legs** for flexibility in accommodating seating surfaces
- **Anatomical cushions** to improve comfort and make the walker easier to steer.
- **Electrically-driven** stepless height adjustment.
- **Leg supports** for extra stability while standing.
- **Rise & Go Harness** provides optimum support for a natural standing action.
- **Ergonomic handles** are adjustable in every direction.
- **330 lb.** weight capacity.

SPECIFICATIONS

Width Legs Open	41.73 inches
Width Legs Closed	26 inches
Weight Capacity	330 lbs.
Length	39.76 inches
Height	35 - 51 inches
Armrests Dimension	13 - 26 inches

EZ56309 - Rise & Go™ w/ Articulating Legs



Rise & Go ambulation Training Harness (Optional Accessory)

The optional ambulation training harness is used for increased safety when training using the Rise & Go. The product's two groin straps and attachment hooks minimize the fall risk if the user becomes unstable while standing or walking (Only for use with the Rise & Go Walker).

"Your Total Patient Lift Solution"

⚠️WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov.

Features



STEP 1

The harness is fitted over the patient's lower back and bottom, the straps adjusted and the leg supports placed in the correct position.



STEP 2

The patient holds onto the handles and the electrically driven power rise function begins. The patient follows a natural pattern of movement with the aid of the harness and the vertical lift action.



STEP 3

On completion of the power rise, the leg supports are swung away and the patient is free to walk. Alternatively, patients place their feet on the optional footplate and are ready to be transferred once standing.

Comfort

Anatomical cushions provide relief for shoulders, arms, and neck.

Electric Adjustment

Electrically-driven height adjustment.

Power Rise

A unique aid for assisting patients and caregivers while standing.

Stirrup for Stability

Stirrup ensures the walker remains stationary while raising the patient.



Adjustable Handles

Ergonomically-designed handles are adjustable in every direction for optimum individual adaptation.

Rise & Go Harness

Provides optimum support, which results in a natural standing action.

Adjustable Leg Supports

The leg supports provide comfort while standing, are height adjustable, and can be moved aside while walking.

Articulating Legs

Flexibility in accommodating seating surfaces