



EZ Way, Inc.

# Limb Strap Instructions

The Limb Lifting Strap is to be used with any EZ Way Lift or EZ Ceiling Lift to provide access for wound care or hygiene or to provide support to a limb that otherwise would have to be lifted and held manually.

1. To use the strap, place the strap under the patient's limb that you would like to raise.
2. Position the strap at the point on the limb where you would like lift. Recommended placement is between the patient's knee and ankle when lifting a leg or elbow and wrist when lifting an arm. NOTE: Do not place the strap below the ankle for lifting the leg.
3. Hook each loop end of the lifting strap to one of the two hooks on the lift hangar bar assembly, hooking one end of the strap on each side. NOTE: Be sure to attach to the center hooks on the hangar bar assembly, NOT the swiveling hangar bar hooks.
4. Use the EZ Lift to raise the limb to the desired height. Use caution when lifting. To ensure the strap does not slip, do not over-extend the patient's limb.
5. Do not leave the patient unattended while using the limb lifting straps.

(Lower limb lift below knee)



(Lower limb lift above ankle)



(Limb lift below elbow)



EZ Way, Inc.  
701 E. Washington St.  
P.O. Box 89  
Clarinda, IA 51632  
[www.ezlifts.com](http://www.ezlifts.com)

1-800-627-8940 or (712) 542-5102  
Fax (712) 542-1899