



EZ Way, Inc.
710 E. Main St.
P.O. Box 89 • Clarinda, IA 51632
800-627-8940 • Fax 712-542-1899
www.ezlifts.com

Instructions for Harness with Leg Straps



1. With the Harness with Leg Straps in hand, lean the resident/patient forward (while you are standing in front of them) and drop the Harness down to approximately 2in below their arms. Securely fasten the safety strap around the patient's chest to a snug fit.
2. With one hand guide and tuck the leg strap around and slightly under their thigh. Then, while raising their leg slightly, tug the leg strap under their upper thigh and rest the strap over their knee. Repeat on the other side.
3. Bring the Stand to the resident/patient.
4. Hook up the harness, using the loops on each side of the harness, to the hooks at the end of the Stand arms. Then hook each leg strap up to the two small hooks located on the neck of the Y-shape that forms the base of the Stand arms. Use whichever loop strap is reachable. You do not have to cross the straps in front.
5. Make sure that the loops are securely placed in the hooks before you engage the motor. Raise the resident/patient.

Questions? Please call 1-800-627-8940