



EZ Way, Inc.

# Hourglass Sling Sizing Chart



## Sling Color Coding System

<b>Silver</b>	<b>Gray</b>	<b>Beige</b>	<b>Burgundy</b>	<b>Green</b>	<b>Black</b>	<b>Brown</b>
All Sizes Antimicrobial	S	M	L	XL	XXL	XXXL

Sling Size	Small	Medium	Large	XL	XXL	XXXL
<b>Weight of Patient</b>	70-100 lbs.	90-220 lbs.	190-320 lbs.	280-450 lbs.	400-600 lbs.	600 + lbs.
<b>1. While log-rolled, maximum distance from back of knee to base of neck</b>	42"	53"	54"	58"	58.2"	60"
<b>2. Maximum upper leg breadth (measured side of right thigh around the buttocks to side of left thigh)</b>	36"	40"	42"	44"	44"	53"
<b>3. Maximum torso breadth (measured around back side of body at the bottom of the rib cage)</b>	24.5"	24.5"	27"	28.5"	37"	42"
<b>4. Maximum Shoulder breadth (measured horizontally between points of maximum protrusion of the deltoid muscles on the upper, outer border of the arm and shoulder)</b>	33"	39.5"	51"	55"	55.5"	55.5"

**NOTE:** The size/weight designations are merely estimates and basic guidelines. A proper fit will depend on factors other than weight measurements, including the height and girth of a patient. A proper fit will involve the judgment of the caregiver. Please consult your EZ Way sales representative for questions.

- It is important to evaluate the width of a patient in relation to the width of the sling, particularly Steps 2 and 4. It is important that no portion of the patient overlap the sides of the sling in these two areas.
- The built in head support should not be below the shoulders.
- **Note:** Slings using Wipeable fabric use beige binding for all sizes and do not follow the color coding system. Slings using Antimicrobial fabric use Silver binding for all sizes and do not follow the color coding system.

