

The Seat Strap is an optional accessory that is used for additional lower body support and can be used for transferring or ambulation. Never use the seat strap without the harness. The seat strap is not a lifting accessory and should not be so taut that it lifts the patient when raising or lowering.

500434  
Seat Strap- Regular



500435  
Seat Strap- 800 lb. Smart Stand



### Applying the Seat Strap

1. Before using the Seat Strap, first fit the harness to the patient, then attach the harness to the stand in the normal fashion.
2. Extend the Seat Strap to its fullest length. Attach one seat strap loop (at the end of the Seat Strap) to the same hook the harness is attached to, located at the end of the stand arms.
3. Make sure the Seat Strap is loosely placed on the backside of the patient. Position the padded side of the seat strap so it faces the patient.
4. Repeat Step 2 with the other seat strap loop on the opposite side, attaching it to the same hook the harness is attached to, located at the end of the stand arms.
5. Have the patient place their arms outside of the straps and place their hands on the handles (or stand arms when using the Classic Stand).
6. Stand beside the patient, and using the hand control, press the UP button to raise the patient slightly off the surface. Tighten harness buckle.
7. Slide the Seat Strap under the patient's buttocks.
8. Press the UP button, and raise the patient to a standing position.
9. Adjust the Seat Strap to the desired tension.

