


## How to Use the Multi-Purpose Sling

Who to use the multi-purpose sling on?

1. Patients who are either single or double amputees above the knee.
2. Someone with very large, fleshy thighs.

How to use:



When using the multipurpose sling, make sure that the labels are away from the patient's body and the center of the commode hole opening "D" is one inch below the Tailbone. After the sling is under the patient, begin hooking the sling to the lift from the head end first. Hook the "A" straps to the lift like a normal EZ sling. Using the different loops on the strap will adjust the way the patient sits on the sling (i.e. upright or reclined). On the thighs there are two ways to hook up the sling. If the amputation is extremely short you must take the right side, thigh loop "C" and go under both legs to the left "B" strap. Feed the "B" strap through the "C" loop, then hook the "B" strap to the lift at the proper height. Repeat this for the other side. The Second way goes between the both legs to achieve a sling look. Bring the right "C" loop under the right thigh, making sure that the sling is smooth., then back over the thigh to the right "B" strap. Feed the right "B" strap through the right "C" loop, then hook up to the lift at the proper height. Repeat for the other side. If you have any questions, please call Sales Support at 1-800-627-8940.

EZ Way, Inc. • (800) 627-8940 • Fax (712) 542-1899 • [www.ezlifts.com](http://www.ezlifts.com)

**⚠WARNING:** Cancer and Reproductive Harm - [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Form #2-006  
08/09/2018