



EZ Way, Inc.

QUICK REFERENCE GUIDE SMART STAND



This quick Reference guide has been created to assist healthcare professionals and other users with only the functional aspects of the system. The guide does not replace the general safety warnings, cautions and recommendations provided in the product instructions for use the **EZ Way Smart Stand**.



REFER TO THE PRODUCT INSTRUCTIONS FOR USE BEFORE USING THIS PRODUCT.

Product Overview



Attach harness loops to “pigtails”.



Use the hand control to operate.





Only for emergency, use the emergency lowering handle to lower the person.

NOTE: The footplate can be removed to ambulate. If using the Smart Stand for ambulation, the rear wheels must be locked when raising or lowering the person. Connect the harness using the long loops.



WARNING: Do not remove the shin pad when ambulating. Removing could result in injury to the user. When refitting, ensure the footplate is reinserted properly.

1. Fit the the harness to the resident, ensuring correct size is used. Ensure the safety strap is buckled and fastened securely around the person’s torso and the person’s arms are outside the harness.
2. Position the Smart Stand in front of the person, open the chassis legs if necessary, have person place their feet on the footplate with their knees touching the shin pad. Do not lock rear wheels when using the footplate.
3. Lower the stand arms to lowest position to aid in hooking up the harness, and adjust shinpad if necessary.
4. Attach the harness loops to the stand arm “pigtail” using the shortest loops when possible.
5. Using the hand control, raise the person to a standing position while tightening the safety strap around the torso while they are raised.
6. Close-up legs and transfer person to desired location.
7. Position the person in front of chair, bed or toilet and lower to a sitting position while standing beside the person.
8. Unhook harness loops from the “pigtails” and move the unit away from the person making sure they lift their feet of footplate.
9. Unbuckle the safety strap and remove the harness from around the resident.



**Training Videos
for Smart Stands**

<https://www.ezlifts.com/ez-smart-stand-training>