



*“Your Total Safe Patient Handling Partner”*

**EZ Way, Inc.**

## EZ Way Accessory Inspection Checklist

The EZ Way slings and harnesses are made specifically for EZ Way lifts and stands. For the safety of the patient and caregiver, only EZ Way accessories should be used with EZ Way lifts and stands. Before operating the unit, ensure the accessory is not ripped, frayed or showing signs of wear.


- 1** Check binding and loops for any fraying, wear, nicks, or tears. Replace if fraying or wear is found.

**EXAMPLES OF ACCESSORIES THAT NEED REPLACEMENT**



- 2** Check sling or harness body fabric for any rips, holes, fraying, or weak spots. Hold sling up to a light. The sling must be replaced if any light shows through.

- 3** Check all stitching of strap connections, (Box-X stitching) for loose stitches. Check all binding stitching for loose stitches. Pull on all straps/loops in opposite direction and note if stitching become loose or comes apart. Remove from service if stitching tears.



- 4** Check handles for loose stitches and tears. Remove from service if any found.



- 5** Check all straps where they attach to the sling hanger bars for wear. Remove from service if any fraying is found.



EZ Way warrants that our slings and harnesses will be free from defects in material and workmanship for a period of six months (see Limited Warranty for full description). Slings and harnesses have a life expectancy of up to two years depending on use, care, and proper laundering process. Slings or harness should be discarded and replaced if showing any sign of damage or wear.

PO Box 89, 710 E. Main. Clarinda, IA 51632 | 1-800-627-8940 | sales@ezlifts.com