



EZ Way, Inc.

## QUICK REFERENCE GUIDE SMART LIFT



This quick Reference guide has been created to assist healthcare professionals and other users with only the functional aspects of the system. The guide does not replace the general safety warnings, cautions and recommendations provided in the product instructions for using the **EZ Way Smart Lift**.



**REFER TO THE PRODUCT INSTRUCTIONS FOR USE BEFORE USING THIS PRODUCT.**

### Product Overview



Attach sling loops to the hanger bar hooks using the same length and color of loop on each side.



Use the hand control to operate.



Only for emergency, use the emergency lowering handle to lower the person on the 500 lb. & 600 lb. models. Turn the red collar clockwise on the 1,000 lb. lift.



1. Fit the the sling to the resident, ensuring correct size is used. Ensure the sling is attached properly following the instructions for that sling.
2. Position the Smart Lift in front of the person, open the lift legs, if necessary, to get around obstacles.
3. Lower the boom and position green nose cone over person's torso. Do not lock the rear wheels.
4. With the hanger bars parallel to the person's body, attach the sling loops to the hanger bar hooks using the same length and color of loop on each side.
5. Using the hand control, raise the person to to the needed height for the transfer.
6. Close-up legs and transfer person to desired location.
7. Position the person over the chair, bed or toilet and lower to a sitting position or supine position, dependent on the transfer.
8. Lower the boom enough to remove the sling loops from the hanger bar hooks and move the unit away from the person.
9. Remove the sling from around the person, unless the sling designed to remain under the person.






Figure 3

- 1T – 1st Loop (shortest)
- 2T – 2nd Loop (middle)
- 3T – 3rd Loop (longest)
- 4 – Center Handle
- 5 – Side Handles
- 6L – Left Sling Leg

- 6R – Right Sling Leg
- 7 – Base of Sling
- 1B – 1st Loop (shortest)
- 2B – 2nd Loop (middle)
- 3B – 3rd Loop (longest)



**Training Videos  
for Smart Lifts**

<https://www.ezlifts.com/ez-smart-lift-training>