

The Chest Strap is used with EZ Way sit-to-stands to help keep forward-leaning people upright and to prevent them from slumping forward. The Chest Strap can be attached to the black padded handles on EZ Smart Stands or to the stand arms on the Classic Stand.

Attaching the Chest Strap to the EZ Stand

1. Lower the stand arms to the lowest position.
2. With the strap-side of the chest strap facing up and the label visible, grasp the shortest strap on either side of the pad.

Steps 3 & 4 for EZ Smart Stands

3. For EZ Smart Stands, approach the lower handles of the EZSmart Stand to attach.
4. With each buckle/receptacle in hand, reach around each lower handle and buckle the short straps behind the chest strap. Tighten the strap accordingly. See FIG 1.

Steps 3 & 4 for Classic Stands

3. For EZ Classic Stands, approach the stand arms to attach.
4. With each buckle/receptacle in hand, reach around the stand arms and buckle the short straps behind the chest strap. Tighten the strap accordingly. See FIG 2.

5. Follow the EZ Smart Stand or Classic Stand operating instructions to attach the EZ Way harness and raise the person in the stand.
6. Once a person is in a standing position using the harness, the caregiver may position and buckle the second set of straps behind the person, over the harness. Use of these straps around a person is an optional procedure.

