



EZ Way, Inc.

"Your Total Safe Patient Handling Solution"

EZ Way Single-Patient Use Accessories

Our single-patient use accessories are an effective tool in helping with infection control. They are designed for use with a single-patient and can be reused on the same patient until they are discharged or the accessory becomes soiled. Available in a variety of styles.



SPS200 - Medium (650 lb. cap.)
SPS300 - Large (650 lb. cap.)



SPSHD100 - Small with Head Support (650 lb. cap.)
SPSHD200 - Medium with Head Support (650 lb. cap.)
SPSHD300 - Large with Head Support (1,000 lb. cap.)
SPSHD400 - XL with Head Support (1,000 lb. cap.)
SPSHD500 - XXL with Head Support (1,000 lb. cap.)



SPH200 - Medium Harness (650 lb. cap.)
SPH300 - Large Harness (650 lb. cap.)
SPH400 - XL Harness (800 lb. cap.)



MULTIPURPOSE
SPSMP300 - Large (650 lb. cap.)
SPSMP400 - XL (1,000 lb. cap.)



SPLS300 - Limb Lifting Strap (300 lb. cap.)



SPTR100 - Turning/Repositioning Sheet (1,000 lb. cap.)

Common Applications:

- Single-patient accessories are 100% breathable and made from a durable yet comfortable polyester material.
- Only for use as a single-patient accessory. Product should not be recycled or redistributed. Single-patient accessories are not designed for long-term use.
- Do not launder. If the accessory is laundered, it should be discarded.
- Operator must inspect accessory prior to each use. The accessory should be replaced at any sign of wear, tear, or soiling.
- Weight load must only be used to hold the amount of weight dictated by the capacity of the lift or stand and the accessory. Equipment and accessory should be used by trained personnel in accordance with operation manual, video, and training checklist to avoid injury.

EZ Way, Inc. • PO Box 89, 807 E. Main • Clarinda, IA 51632 • 1-800-627-8940