



EZ Way, Inc.

"Your Total Safe Patient Handling Solution"

EZ Way Harnesses

Used with EZ Way sit-to-stands, our harnesses are designed for maximum comfort including padded arm wings on our reusable models. Available with or without leg straps, and in a variety of fabric choices. All reusable harnesses have a 1,000 lb. weight capacity.



Regular Harness

- 50945** - Small
- 50910** - Medium
- 50925** - Large
- 50970** - XL
- 50983** - XXL
- 50989** - XXXL



Wipeable Harness

- 500423** - Small
- 500424** - Medium
- 500425** - Large
- 500426** - XL
- 500427** - XXL
- 500428** - XXXL




Antimicrobial Harness

- 500613** - Small
- 500614** - Medium
- 500615** - Large
- 500616** - XL
- 500617** - XXL
- 500618** - XXXL



Anti-Slip Harness

- 50945N** - Small
- 50910N** - Medium
- 50925N** - Large
- 50970N** - XL
- 50983N** - XXL



Single-Patient Use Harness

- SPH200** - Medium Harness (650 lb. cap.)
- SPH300** - Large Harness (650 lb. cap.)
- SPH400** - XL Harness (800 lb. cap.)



Regular Harness with Leg Straps

- 50398** - Small
- 50406** - Medium
- 50408** - Large
- 50363** - XL
- 50385** - XXL
- 50386** - XXXL



Wipeable Harness with Leg Straps

- 500429** - Small
- 500430** - Medium
- 500431** - Large
- 500432** - XL
- 500433** - XXL
- 500473** - XXXL

Fabric Applications:

- Single-patient use harnesses help to reduce the risk of cross-contamination while increasing the level of infection control. 100% breathable and made from a durable yet comfortable polyester material. Only for use as a single-patient accessory. Do not launder or redistribute. Single-patient accessories are not designed for long-term use.
- Our wipeable fabric has been tested and approved with most standard germicidal cleaners.
- Our antimicrobial polyester fabric is made with an antimicrobial treatment which resists bacterial growth, mildew, rot, and provides dust mite protection.
- Our regular harness is made of durable polyester fabric designed to be laundered.
- The neoprene-like fabric of our anti-slip harnesses helps to prevent it from slipping up on the patient.
- Weight load must only be used to hold the amount of weight dictated by the capacity of the stand and the harness. Equipment and harness should be used by trained personnel in accordance with operation manual, video, and training.

EZ Way, Inc. • PO Box 89, 807 E. Main • Clarinda, IA 51632 • 1-800-627-8940

⚠WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov.

Form #2-378
Rev. 11/10/2022