




EZ Way, Inc.

"Your Total Patient
Lift Solution"

Pelican Sit, Slide & Stand Pad

A one way slide pad to help prevent sliding forward in a wheelchair.



Applications:

- Prevents resident from slowly sliding forward in a wheelchair
- Helps to maintain a good sitting position in the wheelchair

Features and Benefits:

- Allows a person to sit down at the front of a chair, and then slide their hips back into the chair
- Webbing handles allow one or two nurses to assist in sliding the resident to the back of the wheelchair
- Special material allows the pad to only slide in one direct and prevent a resident from slowly sliding forward.
- The pad does not build up heat, is comfortable to sit on and is virtually unnoticeable when in use



Item #407H- 15" wide x 22.65" deep

PO Box 89, 807 E. Main • Clarinda, IA 51632 • 1-800-627-8940

Form #2-327
Rev. 06/09/19

⚠WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov.