



EZ Way Hourglass Slings

The hourglass-shaped sling is designed for maximum comfort including special head support. Available with solid or mesh material and 'toileting hole' option.



Common Applications:

- Used with your Patient/Resident where the sling would be left underneath them while sitting in a chair/wheel chair.*
- Amputee above the knee, either single or double.
- Patients with sensitive skin, brittle bones and other special conditions which might apply.
- Patient/Resident wearing a dress or skirt and dignity prevails.
- Not for use with combative or uncooperative Patient/Resident.



Mesh Hourglass Sling w/ Toileting Hole

Form 2-288 Rev. 7/31/18

^{*}The hourglass sling is manufactured in different materials with special options. If this is to be left under a Patient/Resident for a period of time, EZ Way recommends using the mesh hourglass sling. Ask your EZ Way Territory Sales Manager or call 800-627-8940 for more information.