


# **Complete Shower System with Thermostatic Mixing Valve**

"Your Total Patient Lift Solution"



2810 - Shower Panel



The 2810 is a shower panel with thermostatic mixing valve and a separate cleaning system. Pairing it with EZ Way shower trolleys creates a complete shower system that is comfortable, safe and hygienic for both resident and attendant.

#### **Features**

- Easy to operate.
- Intergrated hygienic and safe cleaning system for easy cleaning of equipment.
- Built in bottle for cleaning detergent which are automatically mixed with cold water in the cleaning system.
- Mechanical lock for cleaning detergent bottle.

## Hot and cold water supply for patient shower:

- Thermostatic mixer
- Valve for patient shower
- Shower spray handle
- Patient shower hose white

## **Cold water supply, for cleaning:**

- Safety valve for equipment cleaning.
- Cleaning spray handle. NOTE: Never to be used on patients.
- Cleaning shower hose yellow.
- Cleaning detergent filling.
- Cleaning flow setting valve inside panel.

#### **Options**

• Thermometer on the shower hose to read shower temperature.

Width	23"
Height	15"
Depth	5"
Recommended mounting height from floor	49"
Min. Dynamic water pressure Max. Dynamic water pressure	0.15 Bar 5 Bar