



EZ Way, Inc.

EZ Way Harness Sizing Chart



Harness Color Coding System

Silver	Gray	Beige	Burgundy	Green	Black	Brown
All Sizes Antimicrobial	S	M	L	XL	XXL	XXXL

(Applies to washable and disposable harnesses.)

Harness Size	Small	Medium	Large	XL	XXL	XXXL
Weight of Patient	70-100 lbs.	90-220 lbs.	190-320 lbs.	280-450 lbs.	400-600 lbs.	600+ lbs.
Circumference of patient's torso where harness is applied	26" - 38"	34" - 46"	40" - 56"	50" - 64"	55" - 72"	65" +

- ◆ **Note:** The size/weight of designations are merely estimates and basic guidelines. A proper fit will involve the judgement of the caregiver.
- ◆ Color coding applies to the binding of the harnesses.
- ◆ **Note:** Harnesses and accessories using Wipeable fabric use burgundy binding for all sizes and do not follow the color coding system. Harnesses using Antimicrobial fabric use Silver binding for all sizes and do not follow the color coding system